

Concepts of Imaging and Knobology

Pravin Patil, MD FACC FASE
Associate Professor of Medicine
Director, Cardiovascular Disease Training Program
Lewis Katz School of Medicine at Temple University

Disclosures

- ✓ No relevant financial disclosures



Echocardiography


- ✓ Tomographic Imaging
 - Echo is a “thin slice” imaging tool like cardiac CT, MRI and nuclear imaging
 - However, images are not automatically acquired
 - Achieved by probe manipulation, acoustic windows, patient positioning, balancing of artifacts and image processing
- ✓ The ability to optimize the image acquisition and processing is part of competency in echocardiography

“Knobology”



Source: Philips, GE, Siemens, Esaote, Toshiba


Ultrasound Probe



Source: WikiRadiography.net





Phased Array



Garbi, M. The EAE Textbook of Echocardiography 2011

Machines and Knobs






Rise of the Touchscreen and Trackpad





Source: Philips Healthcare, GE Healthcare

Echo Screen Anatomy




Echo Screen Anatomy



Frame Rate

- ✓ Temporal Resolution
 - Limited by line density and time to scan
 - Improves with narrow sector width
 - Unless line density decreases
- ✓ Appearance
 - FR < 15 Hz appear choppy
 - FR > 15 Hz appear smooth




Echo Screen Anatomy




Imaging Depth

✓ Imaging Depth


- Start with larger than needed
- Adjust to place ROI approximately $\frac{3}{4}$
- Leave small area behind to observe useful artifacts like shadowing




Depth





Depth




Echo Screen Anatomy






Axial Resolution \times Frequency



Zoom



Larger Pixels, Higher Frame Rate
Not Higher Axial Resolution




Harmonic Imaging


- Improves signal-to-noise ratio
- Contrast
 - Non-linear resonance of bubbles to compressions and rarefactions of ultrasound wave
- Tissue (incidental discovery)
 - Related to propagation of sound through the myocardium
 - Non-linear response due to higher speed during compression than in rarefaction.

McCulloch, et al. JASE 2000

Tissue Harmonic Imaging

- ✓ Non-linear distortion of acoustic signal in tissue generates harmonics
- ✓ Noise/artifacts generate no significant harmonic
- ✓ Tissue Harmonic Imaging takes advantage of increased SNR






Fundamental


Second Harmonic

Garbi, M. [The EAE Textbook of Echocardiography](#) 2011

Tissue Harmonic Imaging




Fundamental




Tissue Harmonic




Tissue Harmonic Imaging




Bubbles Have Harmonics too..



Echo Screen Anatomy



Lateral Resolution



- Focal Zone / Focus
- Increased ability to discern two separate objects at the FZ

Apical Wall Motion Abnormality,
Concern for LV thrombus



Echo Screen Anatomy



Gain





Compression





Adapted from Garbi, M. [The EAE Textbook of Echocardiography](#) 2011

Auto Gain





2D Gain




Decrease


Doppler Gain




Time Gain Compensation (TGC)




Lateral Gain Compensation (LGC)



Compression




Color




Human Grayscale

- Humans perceive approximately **30 shades** of gray
- Cones are high bandwidth and allow humans to see at least **10 million shades** of color



Echo Screen Anatomy




Ultrasound Bioeffects

- ✓ Ultrasound is mechanical energy
 - Thermal effects
 - Mechanical effects
- ✓ Background
 - No evidence that diagnostic ultrasound produces harm
 - Subtle or transient effects not well understood
- ✓ ALARA
 - As Low As Reasonably Achievable

Mechanical Index (MI)

- ✓ Quantification of US acoustic intensity
- ✓ $MI = P(\text{Pascals}) / \sqrt{\text{Frequency}(\text{MHz})}$
- ✓ Non-thermal (Mechanical) Bioeffects
 - MI expresses the acoustic pressure of US beam on insonated structures
- ✓ Lower MI induces increased bubble resonance and harmonics

Low MI Example



Thermal Index (TI)

- ✓ Quantification of potential for tissue heating
 - As ultrasound travels through tissue energy is absorbed by tissue and converted to heat
 - Frequency and intensity dependent
- ✓ Recommendation is to keep tissue heating $< 1.5^{\circ}\text{C}$
 - Caution with the febrile patient
- ✓ Thermal Bioeffects

Doppler Echocardiography

- ✓ Optimal 2D images when ultrasound beam is perpendicular to structures
- ✓ Optimal Doppler imaging when ultrasound beam is parallel to flow
- ✓ Apical views allow alignment with most cardiac flows (i.e. aortic, mitral and tricuspid valves)




Doppler Echocardiography

- ✓ Color Flow Doppler
 - Pulse wave modality that cannot resolve high velocities
 - Turbulence/variance maps can help define jet, direction and turbulence
- ✓ Pulse Wave Spectral Doppler
 - Range specific
 - Subject to aliasing at high velocities like CFD
- ✓ Continuous Wave Spectral Doppler
 - Able to resolve high velocities
 - Range ambiguous




Color Flow Doppler


- Pay attention to the baseline
- Make note of the Nyquist limit
- Color scales vary
- Variance maps (see example)
- Optimize size of box and sector for frame rate





Variance Map



Spectral Doppler Scale




Spectral Doppler Baseline





Baseline Adjusted





Doppler Gain




Doppler Gain




Nyquist Limit



Doppler Sweep



75 mm/s



150 mm/s



ASCeXAM Focus

- ✓ Image optimization “Knobology” is part of achieving competency in echocardiography
- ✓ Key components of image optimization
 - 2D Gain
 - Compression
 - Mechanical index
 - Harmonic Imaging
 - Frame Rate
 - Resolution
 - Depth
 - Doppler Gain
 - Doppler Sweep/Scale



